## BACtrack Mobile SMARTPHONE BREATHALYZER

The award-winning BACtrack Mobile integrates a smartphone app and professional-grade breathalyzer to bring blood alcohol content (BAC) wirelessly to your device. Quickly and easily estimate your BAC, monitor your results over time, and with ZeroLine®, estimate when your BAC will return to 0.00%. You can even call an Uber™ directly from the BACtrack App if you've had too much.



INTENDED USE: Screens for the presence of ethanol in human breath.



#### Key Features

- **Xtend® Fuel Cell Sensor** The same professional-grade technology trusted by hospitals, clinics, and law enforcement
- Ultra-Portable Design Mobile weighs just 2oz and easily fits in your pocket or purse
- Bluetooth Connectivity Wirelessly connects to your smartphone via Bluetooth.
   Compatible with most Apple, Samsung, and Google devices
- **ZeroLine® Technology** Estimates when your BAC will return to 0.00%
- BACtrack Companion App Saves and tracks BAC results over time



# BACtrack Mobile SMARTPHONE BREATHALYZER

#### Easily and Accurately Estimate Your Alcohol Level in Seconds

For those seeking the highest level of accuracy and consistency, BACtrack Mobile features Xtend® Electro-Chemical Fuel Cell Sensor Technology and simple one-touch operation. A user simply blows through the mouthpiece and their estimated Blood Alcohol Content (BAC) results are displayed in the app withinin seconds.

### Benefits of Xtend® Fuel Cell Sensor Technology

- Professional-Grade Accuracy
- Precise, Consistent BAC Results
- Long-term Reliability
- Sensitive Enough for High BAC Testing
- Faster Warm Up Time

**BAC**track

- Longer Battery Life
- No False Positives



#### **Product Specifications**

Dimensions
Weight
Sensor Technology
Detection Range
Memory
Power Supply
Warm-up Time
Blowing Time

1.75 x 0.63 x 2.75 in.
2 oz. (57g)
Xtend® Platinum Fuel Cell
0.000-0.400% BAC
Saved in BACtrack App
Mini-USB Rechargeable
10-20 seconds
5 seconds

BACtrack 300 Broadway, Suite 26 San Francisco, CA 94133

